

农业元素分析检测领域最新应用进展-

纳米颗粒

大数据分析

余晶晶

ICP-MS应用工程师

纳米颗粒



纳米颗粒物简介

来源



自然产生


风化，矿化，火山喷发，沙尘暴.....



人类活动

排放到大气、水、土壤中

石化燃烧，机械磨损，金属冶炼和产品生产，机动车.....



人工合成纳米材料 (ENPs)


根据需求而特别设计的纳米材料

碳纳米管，电子开关，药物...

纳米颗粒物简介

人工合成纳米材料 (ENPs)

- 日用消费品-防晒霜，化妆品...
- 食品添加剂和包装
- 药物和靶标
- 消毒纱布
- 半导体生产
- 超导材料
- **农药**
- ...



Many Methods for NP characterization

- **Imaging** methods (TEM, SEM, AFM)
- **Spectroscopic/optical** methods (UV-Vis, dynamic light scattering)
- **Hyphenated techniques** (Chromatographic (or other online) separation coupled with ICP-MS detection).

FFF-ICP-MS

CE -ICP-MS

HPLC-ICP-MS



- **Single particle ICP-MS – today's topic**

- ◆ Particle concentration – mass ($\mu\text{g/L}$)
- ◆ Particle number (particles/L)
- ◆ Particle size and size distribution
- ◆ Particle composition


ICP-MS can perform all these measurements in a single analysis

What is single particle (sp)-ICP-MS


何为单颗粒ICP-MS



常规 ICP-MS: 元素浓度及同位素比分析



单颗粒(sp)-ICP-MS: 纳米颗粒分析



Case : Detection of CuO nanoparticles in edible plants

案例：可食用植物中氧化铜纳米颗粒分析

In collaboration with:




J Nanopart Res (2018) 20: 101
<https://doi.org/10.1007/s11051-018-4192-8>

RESEARCH PAPER

Detection of nanoparticles in edible plant tissues exposed to nano-copper using single-particle ICP-MS

Arturo A. Keller • Yuxiong Huang • Jenny Nelson



- Rinse with DIW
- Spike 20-100 nm CuO NP
- Dry for 2 hours
- Rinse with DIW and collect
- Digest with Macerozyme R-10 (25 °C for 24 h)
- Dilute and analysis with sp-ICP-MS


[A.A. Keller, Y. Huang, and J. Nelson, Detection of nanoparticles in edible plant tissues exposed to nano-copper using single-particle ICP-MS, *J Nanopart Res*, 2018, 20: 101]

7900 ICP-MS parameters

Parameter	Value
RF power (W)	1550
Carrier gas flow rate (L/min)	0.67
Make-up gas flow rate (L/min)	0
Spray chamber temperature (°C)	2
Nebulizer pump (rps)	0.1
Sample depth (mm)	8.0
Oxide ratio (%)	1
Integration time (µs)	100
Acquisition time (s)	90
Mass monitored	⁶³ Cu
ORS ⁴ cell mode	No gas


Quantifying Copper Nanoparticles
on Plant Leaves using Single-Particle
ICP-MSInvestigating the adsorption and release of nano-
pesticides on plants using the Agilent 7900 ICP-MS

Control samples – leaves not exposed to CuO NPs




A green leaf lettuce, **B** collard green, **C** kale
Surface CuO NPs size distribution

Exposed samples – leaves exposed to CuO NPs



Surface CuO NP concentration




Surface CuO NP concentration

数据界面

-----纳米颗粒分析结果

- # of Particles
- Conc.(particles/L)
- Ionic Conc.(ng/L)
- Background Equivalent Diameter (nm)
- Particle size (nm)

		197 Au							
Type	Sample Name	Nebulization Efficiency	# of Particles	Conc. (particles/L)	Conc. (ng/L)	Ionic Conc. (ppb)	BED (nm)	Particle Size (nm)	
Sample	EtOH10%								
RM	Au 60nm 100ppt	0.069	1377	5.6E+7	100.0	0.0469	5.65	55	
Sample	Au 30nm 10ppt	0.069	988	4.0E+7	9.5	0.0133	3.77	27	
Sample	Au 60nm 100ppt	0.069	1347	5.5E+7	103.2	0.0466	5.73	56	
Sample	60nm 50ppt 30nm 5ppt	0.069	1116	4.6E+7	55.9	0.0362	5.26	53	
Sample	60nm 80ppt 30nm 2ppt	0.069	1227	5.0E+7	84.3	0.0404	5.46	56	
Sample	60nm 20ppt 30nm 8ppt	0.069	755	3.1E+7	25.3	0.0192	4.26	28	
Sample	60nm 50ppt 30nm 5ppt	0.069	1115	4.6E+7	58.0	0.0337	5.14	54	
Sample	60nm 80ppt 30nm 2ppt	0.069	1217	5.0E+7	84.0	0.0433	5.59	56	



Au 60nm 100ppt 197 Au 21 Samples (21 total)

强大的软件功能及数据处理能力

Batch - FTPA-8900s_SUS-02_May_2017-09_12.18b

Save Batch Add to Queue Validate Method Set Q1/Q2 Masses Select RM Mass Tune Mode: <All> Autosampler Nebulizer Pump Speed

Acq Method Data Analysis Method Sample List

Acq Parameters PeriPump/ISIS Tune

Acq Mode

Single Particle Analysis

Tune Mode	#1: Warm NH3 Fe	#2: Warm NH3 Cr	#3: Warm NH3 Ni	#4: Warm NH3 Mn	#5: Warm NH3 Mo
Stabilization Time [sec]	0	5	5	5	5
Scan Type	MS/MS	MS/MS	MS/MS	MS/MS	MS/MS

Acq Option

Auto/Semi Auto Tune before Batch

Generate Tune Report

P/A Factor Adjustment

[Advanced Configuration](#)

Total Acq Time

248.000 sec

Element Name	Monitor	Q1 -> Q2	IntegTime /Mass [sec]	Monitor	Q1 -> Q2	IntegTime /Mass [sec]	Monitor	Q1 -> Q2	IntegTime /Mass [sec]	Monitor	Q1 -> Q2	IntegTime /Mass [sec]	Monitor	Q1 -> Q2	IntegTime /Mass [sec]
Si	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A
P	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A
S	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A
Cr	<input type="checkbox"/>	N/A	N/A	<input checked="" type="checkbox"/>	52 -> 52	0.0001	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A
Mn	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input checked="" type="checkbox"/>	55 -> 55	0.0001	<input type="checkbox"/>	N/A	N/A
Fe	<input checked="" type="checkbox"/>	56 -> 56	0.0001	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A
Ni	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input checked="" type="checkbox"/>	60 -> 60	0.0001	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A
Mo	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	N/A	N/A	<input checked="" type="checkbox"/>	95 -> 95	0.0001	<input type="checkbox"/>	N/A	N/A

²⁶
Fe
Iron
55.847

²⁴
Cr
Chromium
51.9961

²⁸
Ni
Nickel
58.6934

²⁵
Mn
Manganese
54.93805

⁴²
Mo
Molybdenum
95.94

Single run collects multi-element sNP data and combines it into a single data batch table



- Save time
- Reduce the number of autosampler movements and sample contamination risk

大数据分析



Agilent ICP-MS 应用于污染物示踪的工具

ICP-MS 定量分析与 MPP 同时使用



Food Profiling Can-Do in Food


轮廓分析在食品中的潜在应用

- 品种鉴定
- 产地判别
- 真伪鉴别
- 年份及采收季节判别
- 加工工艺质量控制
-



为什么用无机元素组分进行产地鉴定？

农产品吸收土壤和水中的金属，进而在其体内富集




由于不同地区土壤和水中微量元素含量不同，因此农产品中痕量金属的富集也会有差异




(一) 应用ICP-MS判定中药丹参产地

各地区丹参样品的Pb同位素比值的比较



应用ICP-MS解析福建省水体和土质对植物的影响

同一地区（中江）采样的丹参、土壤和水样中的Pb同位素比的比较



MPP软件应用--橄榄油产地判定

MPP: Mass Profiler Professional

Sample preparation

- 3g of Olive oil was weighed into 50 mL conical tubes, and diluted to 10g with a Premisolv diluent
- All samples were put in a mechanical shaker for 20 minutes
- Internal standard was “t-ed” in before the nebulizer (Sc and Y were used)



Country origin	n
USA (California)	2
Greece	14
Italy	6
Spain	7

50
elements


Agilent ICP-MS System

**Agilent 单一He模
式测定50种元素；**

**良好的消干扰技术
避免了假阳性结果**



Filter function in MPP is very useful!



clear

separation

The filters in MPP reduces the data analysis labor by extracting the elements with significant differences

中国蜂蜜种类鉴定 文献

----- 中国检验检疫科学研究院， 庞国芳院士实验室

室

JOURNAL OF
AGRICULTURAL AND
FOOD CHEMISTRY

Article

pubs.acs.org/JAFC

Chemometric Determination of the Botanical Origin for Chinese honeys on the Basis of Mineral Elements Determined by ICP-MS

Hui Chen,^{†,§} Chunlin Fan,^{*,†} Qiaoying Chang,[†] Guofang Pang,^{*,†,§} Xueyan Hu,[†] Meiling Lu,[#] and Wenwen Wang[#]

[†]Chinese Academy of Inspection and Quarantine, No. 3 Gaobeidian North Road, Chaoyang District, Beijing 100123, China

[§]College of Environmental and Chemical Engineering, Yanshan University, Qinhuangdao 066004, China

[#]Agilent Technologies (China) Company, Ltd., No. 3 Wang Jing Bei Lu, Chaoyang District, Beijing 100102, China

ABSTRACT: In this work, the potential of mineral elements and chemometric methods as a tool to classify Chinese honeys according to their botanical origin was examined. Twelve mineral elements (Na²³, Mg²⁴, P³¹, K³⁹, Ca⁴³, Mn⁵⁵, Fe⁵⁶, Cu⁶³, Zn⁶⁶, Rb⁸⁵, Sr⁸⁸, and Ba¹³⁷) of 163 Chinese honey samples, including linden, vitex, rape, and acacia, collected from Heilongjiang, Beijing, Hebei, and Shaanxi, China, in 2013 were determined by the ICP-MS method. Principal component analysis (PCA) reduced 10 variables to four principal components and could explain 93.06% of the total variance. Partial least-squares discriminant analysis (PLS-DA) and back-propagation artificial neural network (BP-ANN) were explored to construct a classification model. By PLS-DA, the total correct classification rates for model training and cross-validation were 90.9 and 88.4%, respectively. By BP-ANN, the total correct classification rates for model training and cross-validation were 100 and 92.6%, respectively. The performance of BP-ANN was better than that of PLS-DA. The validation of the developed BP-ANN model was tested by the independent test set of 42 honey samples. Linden, vitex, and rape honey samples were predicted with an accuracy of 100%, whereas one acacia honey was predicted as rape honey with an accuracy of 92.3%. It is concluded that the profiles of mineral elements by ICP-MS with chemometric methods could be a potential and powerful tool for the classification of Chinese honey samples from different botanical origins.

KEYWORDS: honey, mineral elements, ICP-MS, botanical origin, PCA, PLS-DA, BP-ANN



采用 ICP-MS 和 Mass Profiler Professional 软件测定储存温度和包装类型对葡萄酒中痕量金属成分的影响

应用简报

食品安全

作者


Helene Hopfer 和 Susan E. Ebeler
葡萄栽培和酿酒系
加州大学戴维斯分校
美国

Jenny Nelson
安捷伦科技公司
5301 Stevens Creek Blvd.
加利福尼亚州圣克拉拉市
美国

摘要

本实验建立了一种采用 ICP-MS 测定葡萄酒中 20 多种痕量金属的方法，研究了包装以及储存条件对痕量金属成分的影响。实验发现 5 种痕量金属的检测限 (LOD) 介于 0.001 至 0.044 $\mu\text{g/L}$ 之间，定量限 (LOQ) 介于 0.003 至 0.057 $\mu\text{g/L}$ 之间，具体数据随包装和储存条件的不同会发生显著变化。葡萄酒中元素的含量随这两个变量而变化，其中包装类型（盒装式或带软木塞/螺旋盖的瓶装式）对金属成分的影响最大。

MPP 可用于Agilent 所有质谱



Thank you!

关注售前官方微信号
安捷伦视界
最新方案助您一臂之力
丰富互动，精彩纷呈



访问**中国市场热点页面**
获取行业洞见, 聚焦中
国方案
www.agilent.com/chem/china-topic



服务热线：
400 820 3278
800 820 3278